



Argonne
NATIONAL
LABORATORY

... for a brighter future

Pump-Probe EXAFS or Probing Structural and Electronic Excitations with X-rays

Lin X. Chen (ANL)

David M. Tiede (ANL)

Bernhard Adams (ANL)

Stephe Ross (ANL)

Ray Conley (BNL)

Klaus Attenkofer (ANL)

Michael Goerlich

Guvenc Akgul

Stefan Groissmeier

Guy Jennings

Chuck Kurtz

Xiaoyi Zhang

Markus Winterer (Uni. Duisburg-Essen)

Martin M. Nielsen (CMM/Copenhagen)

Thomas Nielsen (CMM/Copenhagen)

Ben Gilbert (LBNL/Berkeley)

Thomas Wong (IIT)

Jean-Francois Genat (UofC)

Henry Frisch (UofC)



U.S. Department
of Energy



A U.S. Department of Energy laboratory
managed by The University of Chicago

Outline

- What kind of information your synchrotron technique may provide
 - Dynamics of glasses
 - X-ray probes structure and electronic state
- What are the basic principles of your technique
 - The excited state: Intermediates, ensembles, & more
 - The pump-probe experiment
 - Reversible/irreversible problems
 - The time resolution
- What are the limitations of your technique
 - Poor Signal/noise ratio
 - Increase detection efficiency/available flux
 - Sample damage
- What kind of sample does one need


Motivation

Dynamics of Glasses

- Reversible- non-reversible processes
- Typical Problems:
 - Melting
 - Phase transitions (RW-DVD)
 - Dynamic behavior (like.....)
 - Charge dynamics in glasses
 - Chemical reactions
- X-ray will probe ?

Motivation: X-ray Probes Structure and Electronic State

Electronic Structure




Absorber Atom

NEXAFS

- Local Electron Density
- Symmetry
- Atomic selectivity

Local Structure




Absorber Atom

EXAFS

- Distance
- Symmetry
- Atomic selectivity
- Number of Neighbors




Global Structure



WAXS/SAXS

- Global Electron Density
- Global Shape
- Movie with Limited Spatial Resolution



Motivation: Electronic Structure, Symmetry and Pre-Edge



T_1 excited state has (d_{z^2} , $d_{x^2-y^2}$) electronic configuration with singly occupied $3d_{x^2-y^2}$ and $3d_{z^2}$ orbitals .
 The energy gap between $3d_{x^2-y^2}$ and $3d_{z^2}$ in the final excited state is ~ 2.2 eV .
 The energy of $3d_{x^2-y^2}$ orbitals shifts up in the final excited state.


Motivation: Electronic States & White Line

- Change of E_{Fermi} : derivative like shape
- Change of density of state: various peaks may appear
- Energy level of empty states (with restrictions)
- Local symmetry changes




Motivation: Local Structure & EXAFS

Time evolution of As-Se glasses after laser excitation (Se-edge):




- “Movie” of local structure
- Atomic-selectivity
- Melting-effects
- Smaller systematic errors (difference measurements!)

Motivation: Global Structure & GI-SAXS/WAXS





- Crystallization (Gi-WAXS)
- Melting (Gi-WAXS)
- Phase transitions (Gi-WAXS)
- Grain growth (for example catalysis)
- Today: first experiments on LB-films


Hexagonal Cylinder



Hexagonally perforated layer




Gyroid : Cubic($Ia3d$)




Basic Principles: The Excited State

- Experiment on ensemble
- Initial excitation coherent (very fast)
- Thermalization to intermediate (longer lifetime; usually probed)
- Relaxation process includes many transient states/intermediates
- De-phasing of ensemble with time




Basic Principles: The Pump-Probe Experiment



- Pump creates Excite state (intermediate) population
- Probe pulse is delayed and probes a mixture of different excited states & ground state
- Multiple probe-pulses can be used to probe the development of a system (movie)
- Signal is averaged over many excitation pulses (ensemble/reproducibility)

Basic Principles: The excitation Mechanism




Selectivity:

- Polarization
- Resonator /Symmetry
- Wavelength
- Combinations

The spectrum (often Lasers)

- Electromagnetic pulses: magnetic/electric excitations (collective excitation)
- THz excitations: vibration/collective modes
- IR-pulses: vibrations/heating/melting
- Optical & UV: bond-braking
- Ultrafast pulses: optical/acoustic phonons

Basic Principles: The Time-Resolution (stroboscopic)





Time-resolution is given by **convolution of excitation and probe-pulse** (stroboscopic experiment) or by **convolution of probe-pulse and detector resolution**

Typical excitation and probe pulses


- Laser excitation: 10fs-10ns
- Synchrotron pulse: 30ps-1ns (APS typically 80ps)
- XFEL: 1fs-200fs
- Plasma sources: 50fs-500fs
- Thermal sources: 1ps-20ps
- Slicing sources: 200fs-5ps

Basic Principles: The Time-Resolution (time resolution of detector)




- Probe pulse covers the full time evolution of interest
- Time spectrum will be “statistically” acquired (less systematic errors)
- Time resolution is independent from timing-jitter between excitation and probe
- Available detection systems:
 - Streak cameras (~1ps, low quantum efficiency)
 - PIN-diodes (~1ps, high quantum efficiency, small solid angle first demonstrations)
 - PIN-area detectors (~1ps, ASIC design in progress)

Basic Principles: The Time-Resolution (time resolution of detector)



- First tests show resolution of about 3ps
- Only available for single element
- Work on multielement system in progress (first 10-20 element system available in 18 month)
- ASIC design in collaboration with HEP (University of Chicago)
- Detector Unit: InGaAs-PIN diodes with 40 μ m active area

Basic Principles: The Experimental Setup



- Multi-Photon counting
- Digital processing
- Optimized solid angle (20%)
- Lytle-detection principle (Z-1 filter, soller-slit, non energy resolving detector)

The Limitations: Sample Damage and Irreversible Processes

Typical damage:

**Sample moving /
exchanging**

- Evaporation
- Oxidation / Reduction
- Degradation


Refurbishing

Coating

**Vacuum
/inert gas**



Reactor:




**SYSTEMATIC ERRORS AND
LARGE SAMPLE QUANTITIES**

The Limitations: Hugely Different Absorption Coefficient for Pump & Probe

Excitation absorption length:
Probe absorption length:


typically 100nm or shorter
typically 10 μ m or longer



POOR SIGNAL / BACKGROUND PROPERTIES

Experimental ways out:

- Gracing incidence Geometry for x-rays normal (Brewster) for excitation
- Dilute sample system: dispersions of nano-particles
- Excitation wavelength tuning



The Limitations: Largely Different Repetition Rates for Pump & Probe


Laser			X-ray (Synchrotron)
Rep-Rate	Power/pulse I _{se}	Spectral range	
[Hz]	[mJ/pulse]		6.5 MHz (APS standard mode)
1-10	50-100	IR/optical/UV	
1K-10K	1-0.5	IR/optical/UV	
100K-10M	0.05-0.01	IR	



**POOR SIGNAL / BACKGROUND
PROPERTIES**

The Limitations & Solutions High Flux Beamlines



Example of two inline undulator beamline (11-ID-D/APS):



- Long straight section
- Multiple specialized (inline) undulators (restriction of energy range)
- Large acceptance angle of beamline optics
- High heat-load optics & front-end
- Variable bandwidth monochromators


Limitations & Solutions: Micro/Nano-Focus Provides Flexibility for Excitation Wavelength

New approach using MEMS-technology?



- Wavelength flexibility costs pulse power
- -> low excitation rates
- Strongly focusing to keep power-density constant

- Variable focus size
- Scannable focus spot
- Large collection efficiency
- Very fast adjustable (shot-to-shot)?
- Beam-shaping capabilities



Limitations & Solutions


The Fluorescence Detection System

APS Mode	24-Bunch	324-Bunch	1296-Bunch
Repetition Rate	6.5MHz	88MHz	352MHz
Bunches/ μ s	6-7	88	352
Photons/Bunch (ID(5×10^{13})/BM(10^{12}))	7.7M / 154K	0.57M / 11.4K	0.14M / 2.8K
Fluo-Photons/Bunch (1%Fe in light matrix)	385K / 7.7K	28.5K / 570	7K / 140
Photons in Detector / Bunch (10% solid angle/ 10% efficiency)	3.8K / 80	285 / 6	70 / 1

A 10.000 element system will be required
to utilize all photons

Limitations & Solutions

The Fluorescence Detection System: Conventional SDD System





- Monolithic front chip (commercially available)
- Preamplifier (linear amplifier) “home-made” and optimized for ADC
- ADC/sender/receiver: industrial collaboration
- Digital processing in consumer electronics (GPU’s?)
- Required time spacing: ~120ns

Cost per channel \$500-\$1000

Total cost about \$1M for 1000-2000 SDD system



Limitations & Solutions

The Fluorescence Detection System: Multi-Photon Detection




- Energy resolution by filter-systems and filters/crystal optics
- Utilization of modern rapid prototyping techniques
- Fast detection systems like plastic-scintillators/PMT or APD's

Filter System with Fe-Fluorescence




Limitations & Solutions

The Fluorescence Detection System: Digital Processing



- Usage of signal averager (1GHz sampling rate / 1ms storage)
- Acts like oscilloscope with very small dead time between triggers (50 μ s)
- Maximal counting time per point about 60s
- Read-out overhead \sim 100ms
- Data analysis happens after memory is read out
- Fitting of the detector response function to the measured signal
- About 5 order of dynamic range
- Maximal 100counts per detector unit (1 photon \sim 25mV peak)



Sample Requirements

- No standard Experiments: Each system has to be discussed
 - Flat thin samples (thickness optimized to optical absorption)
 - Solution/particle systems are possible
- Scientific problem has to be well characterized
 - Optical pre-characterization
 - Large parameter space to probe
 - X-rays are very selective !
- Good knowledge about ground state experiments
 - Sample damage
 - Good reference (for proposed excitation models)
 - Good theoretical knowledge

Conclusion

- Laser initiated time resolved techniques are relative new (since 5 years available)
- Signal/noise ratio & sample damage are the most important limitations
- Many technologies are around to overcome the experimental problems but they are not commercially available (large development effort)
- Typical time resolution is about 80ps (groups are working to reach 1ps)
- X-ray techniques are highly specific probes therefore the scientific problem has to be well defined
- You will need well characterized samples